

Aero

ANTI-DECUBITUS MATTRESS



Features :

This lattice style air mattress is specially designed for bed sore cure. With alternating air pressure, it can effectively curb the growth of microbes. It is economical in cost, simple in structure, easy in operation. It's a perfect instrument made for bed ridden patients suffering from paralysis, burn, after operation etc.

Specification :

Voltage	: AC 220V, 50Hz
Power consumption	: 10 W
Pressure range	: 50 -105 mm Hg
Cycle time	: 12 min
Material	: PVC, 0.3 mm thick
Weight capacity	: 114 KGS (250 lbs)
Rate	: 7 L/Min
Pump size	: 28 x 13 x 10 cm
Pump weight	: 1.565 kgs
Tube length	: 167cm; 2mm thick, 10mm diameter
Tube weight	: 0.31 kg
Mattress size	: 230 x 103 cm (L x W) (before inflation, including flap); 190 x 90 x 6cm (L x W x H) (after inflation, not including flap)
Each flap length	: 48 cm (Each mattress has two flaps.)
Mattress weight	: 2.535 kgs

Imported & Marketed By :



HEMANT SURGICAL INDUSTRIES LTD.

647, Girgaum Road, Avijah Building, Shop No. 6 & 7, Near Metro Cinema, Mumbai-400 002.

Tel. : 2203 0935, 2203 0936 · Fax : 91-22-2207 7585

Website : www.hemantsurgical.com · E-mail : info@hemantsurgical.com