

## AERO DIGITAL THERMOMETER



### FEATURES

Get exact measurement values of the fever with this Digital Classic digital thermometer from AERO

- **Great Accuracy and Range**  
Designed to provide the accurate body temperature, this thermometer has a range of 90 to 109.9 degree Fahrenheit.
- **LCD Display**  
The LCD display is clear and easy to read and shows the accurate body temperature.
- **Memory Function**  
This thermometer stores the value of the last measurement, helping you monitor the trend of the fever.
- **Low Power Consumption**  
The thermometer consumes less power and has approximately 200 hours of battery life.