



ONE MACHINE SERVES TWO PURPOSES

93% High O₂ Purity


Nebulizer function

93% high
O₂ PurityHigh quality
molecular sieveAdjustable
flowrate 1L~9.9L

LED display

Nebulizer
function48 hours continuous
oxygen supplyO₂ Outlet

Nebulizer Port

O₂ Purity %

Clear LED display easy to read



558mm

345mm



280mm

PRODUCT PARAMETER

Model No.	KJR-Y51W
Oxygen purity	93%±3%
Max recommended flowrate	5LPM
Product size	345mm*280mm*558mm
Net Weight	15.5kg
Noise	≤45dB(A)
Outlet pressure	20kPa~70kPa
Power supply	220V±22V, 50HZ±1HZ
Rated power	≤400W