

# INTRAVENOUS ADMINISTRATION

## PEDIATRIC INFUSION SET



With 23G Needle

### Product Features and Specifications

- Non-Vented (Type 200) and Vented (Type 500) are available
- ISO 8536-4 Compliant Spike
  - > Conical design with a graduated shoulder structure allowing for easy closure piercing
- ISO 8536-4 compliant drip chamber (60 drops/ml)
  - > Soft drip chamber facilitates fluid level adjustment
  - > Incorporates a precise dripping port made of stainless steel that provides accurate infusion rate at 60 drops per ml
- Large Size V-Clamp
  - > For better flow control adjustment for pediatric administration
- ISO 594/1 & ISO 594/2 compliant end-connection
  - > Available in luer slip, fixed or free luer lock with cap
- Injection site (Y-injection site, flashbulb or rubber tube)
  - > Isoprene (latex free) injection port
- Transparent Hub
  - > Easy and quick observation of flash back
- Soft Clear Flexible Tubing
  - > Medical Grade PVC, kink resistant for smooth fluid transportation

- 23G Vein Needle
  - > Extra-Thin wall needle for better flow rate with minimal patient discomfort
- Sterile Pack-ETO sterilization according to EN 550 and ISO 11135 standards.

### **Benefits of JMS Pediatric Drip Tip**

- Precise dripping port with stainless steel tip is built into the drip chamber for Pediatric Type Infusion Set
- JMS Pediatric Drip Tip creates smaller drops for accurate control 60 Drops per 1 ml
- Available in Type-200 (Non-vented) and Type-500 (Vented)



**JMS Pediatric Drip Tip**