


TRANSFUSION THERAPY

SCALP VEIN SET - TYPE Y



Scalp Vein Set

Product Features and Specifications

- Super-thin-wall needle provides larger inner diameter that enables sufficient flow and smooth infusion
- Perfect geometry of needle profile that is individually siliconized for sharpness and painless venipuncture
- Unique wing design for easy gripping and fixing on patient arm. Non-slip design allows secure and safe venipuncture
- Optimum softness of luer connectors offers firm connection with luer adaptor, three-way stopcock, syringes..etc
- Luer Lock connector and intermittent configuration available
- Soft Clear Flexible Tubing
 - > Medical Grade PVC, kink resistant for smooth flow and reliable connection
- Individually blister/tape packed using medical grade paper and plastic film
- Compliance to ISO 6009 standards using color codes for needle gauge size identification
- Sterile Pack - ETO sterilization according to EN 550 and ISO 11135 standards
- Soft Clear Flexible Tubing
 - > Medical Grade PVC, kink resistant for smooth transfusion therapy (150 to 180cm)
- 18G Vein Needle
 - > Super-thin wall needle for better flow rate with minimal patient discomfort
- Sterile Pack - ETO sterilization according to EN 550 and ISO 11135 standards

STANDARD

GAUGE (O.D)	NEEDLE LENGTH	WING COLOR	TUBE LENGTH
18G (1.25mm)	3/4" (19mm)	Pink	30 cm
19G (1.08mm)	3/4" (19mm)	Cream	
20G (0.90mm)	3/4" (19mm)	Yellow	
21G (0.81mm)	3/4" (19mm)	Green	
22G (0.71mm)	3/4" (19mm)	Black	
23G (0.64mm)	3/4" (19mm)	Blue	
23G (0.64mm)	5/8" (16mm)	Blue	
24G (0.56mm)	3/4" (19mm)	Red	
24G (0.56mm)	5/8" (16mm)	Red	
25G (0.51mm)	3/4" (19mm)	Orange	
25G (0.51mm)	5/8" (16mm)	Orange	
27G (0.41mm)	1/2" (13mm)	Grey	

INTERMITTENT

GAUGE (O.D)	NEEDLE LENGTH	WING COLOR	TUBE LENGTH
19G (1.08mm)	3/4" (19mm)	Cream	7 cm
21G (0.81mm)	3/4" (19mm)	Green	
23G (0.64mm)	3/4" (19mm)	Blue	
25G (0.51mm)	3/4" (19mm)	Orange	

