

New Latest Technology

Joylife



Air Angel

PEDIATRIC MESH NEBULIZER

Make children enjoy nebulization

BREATH

easy

Hidden medication cup


Eliminate the fear of nebulization for children

Cute cartoon + Breathing light

Attract children's attention

Nearly silent operation, multi-angle nebulization

Suitable for using when children are sleeping.



Small and Uniform particles

Can be absorbed more easily for children



Musical chip built inside

Compatible for downloading children preferred music from computer, making it a pleasant time when use.



Easy to Carry
with convenient carry pouch

ACCESSORIES :

	USB Cable
	Mask
	Mouth Piece

Technical parameters

Model: Air Angel

Power Supply: Internal power supply

Power Consumption: <3.0w

Nebulization Rate: >0.3ml/min

MMAD: 2.5±30%µm

Medication Cup Capacity: 6ml

Device Dimensions: 53(L) * 52(W) * 140(H)mm

Device Weight: 107g



Manufactured and Marketed by :

HEMANT SURGICAL INDUSTRIES LTD.

502, 6th Floor, ECSTASY, City of Joy Commercial Complex, J.S.D. Road Mulund (W), Mumbai - 400 080.

Tel. : 022 - 2203 0936, 2591 2747/5289 • E-mail : info@hemantsurgical.com, sales@hemantsurgical.com

Website : www.hemantsurgical.com