

Neitz Binocular TTL Loupes

BTL-2S / BTL-2N



NEITZ

Made by Neitz for Your Medical Field

NEITZ TTL Loupes

Wide and Clear Field of View with Durability and Safety

● Easy-to-wear TTL loupes with flexible frame.

● Protective shield with scratch-resistant coating.



● Bright, wide and clear field of view. Furthermore, the focus is adjustable.

● For patients' safety, the lens units are fixed with screws to prevent lens unit from dropping.

Specification

Magnification	2.0x to 2.5x
Working Distance (mm)	550 to 350
Field of View Diameter (mm)	113 to 58
Pupil Distance (mm)	BTL-2N: 56 BTL-2S: 62
Declination Angle	40°
Weight (g)	95

*Optional headlight is available. For the detail, please visit our website.



Neitz LED Headlight NSI-X

Neitz Instruments Co., Ltd.

4F Ichibancho Court, 15-21, Ichibancho, Chiyoda-ku, Tokyo 102-0082, Japan
Tel : (+81) 3-3237-0552 Fax: (+81) 3-3237-0554 URL: <https://www.neitz.co.jp/en/>



Neitz website