

ZONCARE-Q7

Full Digital Color Doppler Ultrasonic
Diagnostic System



HEMANT SURGICAL INDUSTRIES LTD

502, 5th Floor, Ecstasy Business Park, Citi Of Joy, JSD Road,
Mulund West, Mumbai, Maharashtra 400080

Email - monil.g@hemantsurgical.com , Marketing@hemantsurgical.c

Contact - +91 8104729765 / +91 9619484952 / +91 7223886262

Website - www.hemantsurgical.com



ZONCARE-07

Full Digital Color Doppler Ultrasonic Diagnostic System

- 19 inch high resolution medical LED monitor
- 10.4 inch Anti-dazzle LED touch screen
- 4 transducer connectors
- Rotated console and foldable monitor
- Proposed USB ports
- Hidden draw-out keyboard
- Friendly transducer bracket design
- Support multiple peripheral interfaces
- Standard configured gel heating cup
- Optional built-in lithium battery for continuously working



Advanced Imaging Technology


Spatial Compounding Imaging
Trapezoidal Extended Imaging
Linear Independent Deflection Imaging
Speckle Noise Removal Technology
Tissue Harmonics Imaging Technology
3D/4D Imaging Technology
Intelligent Navigation Technology
One Touch Optimize Technology
Cartoid IMT Measurement
Dual/Triple Synchronous

Abundant clinical application

Zoncare-Q7 has excellent performance in obstetrics, gynecology, pediatrics, breast, musculoskeletal, heart and blood vessel areas. Build-in perfect analysis package, including heart, obstetrics, urology, abdomen, gynecology, small parts, blood vessel etc.



High Density and Broadband Frequency Transducers



Convex

Linear

Transvaginal

4D volume



Phased

Micro-convex

Puncture guide


Foot Switch

Innovative Ergonomic Design



Advanced Imaging Technology


Speckle Noise Removal Technology




Advanced 3D Imaging



Real-time Triplex Synchronous Imaging



4D Imaging



Pulse Inverse Phase Harmonic Imaging



Trapezoidal Extended Imaging

