



KAI MEDICAL SURGICAL BLADE & SCALPEL



Product specifications may change without notice.


KAI INDUSTRIES CO., LTD.
1110 Oyana, Seki-shi, Gifu 501-3992, Japan
Phone: +81 (0)575 28 5040
Fax: +81 (0)575 28 6611
www.kaimedical.jp/



JAPAN QUALITY
Pursue the safety and quality of life in healthcare globally

DISPOSABLE SURGICAL BLADES

Blade Shape



Specification and Feature

- No.10/No.11/No.12/No.12B/No.15/No.15C fit handle No.3
- No.20/No.21/No.22/No.23/No.24/No.25/Stitch Cutter fit handle No.4
- Blade : Stainless steel / Carbon steel
- Sterilized
- Single use
- 100 pcs/box

Stainless steel



Carbon steel



Size Lineup

Blade Type	Ref No.	
	Stainless steel	Carbon steel
No.10	310-B	210-B
No.11	311-B	211-B
No.12	312-B	212-B
No.12B	312B-B	212B-B
No.15	315-B	215-B
No.15C	315C-B	215C-B
No.20	320-B	220-B
No.21	321-B	221-B
No.22	322-B	222-B
No.23	323-B	223-B
No.24	324-B	224-B
No.25	325-B	225-B
Stitch Cutter	—	2SC-B

DISPOSABLE SURGICAL SCALPELS

Small: No.10/No.11/No.12/No.15



Large: No.20/No.21/No.22/No.23/No.24/No.25



Specification and Feature

- Two types of handle sizes: Large and Small
Small: No.10/No.11/No.12/No.15
Large: No.20/No.21/No.22/No.23/No.24/No.25
- Blade : Stainless steel
- With a transparent cap
- Handle with a scale
- Sterilized
- Single use
- 20 pcs/box

Size Lineup

Blade Type	Ref No.
No.10	510-A
No.11	511-A
No.12	512-A
No.15	515-A
No.20	520-A
No.21	521-A
No.22	522-A
No.23	523-A
No.24	524-A
No.25	525-A



DISPOSABLE SAFETY SCALPELS

Small: No.10/No.11/No.15



Large: No.12/No.20/No.21/No.22/No.23/No.24



Specification and Feature

- The blade can be retracted into the handle.
- The blade extends and retracts with one hand.
- Two types of handle sizes: Large and Small
Small: No.10/No.11/No.15
Large: No.12/No.20/No.21/No.22/No.23/No.24
- Blade : Stainless steel
- Handle with a scale
- Sterilized
- Single use
- 10 pcs/box

Size Lineup

Blade Type	Ref No.
No.10	SS-10
No.11	SS-11
No.12	SS-12
No.15	SS-15
No.20	SS-20
No.21	SS-21
No.22	SS-22
No.23	SS-23
No.24	SS-24



Before Use



After Use

